




Precision Coping Company

TO GET THE BEST RESULTS...



- Ditch your dies directly below the margins
- Don't mark the margins
- Don't use spacer or hardener
- Light, matt die stones scan best
- Send us an opposing model, or a silicone bite taken between the prep and opposing teeth
- When possible, and with complex cases, send us a removable wax-up and we'll follow its contours

30 Clarendon Road
Watford WD17 1JJ

01923 237484

www.precisioncoping.co.uk
info@precisioncoping.co.uk

Precision Coping Company

30 Clarendon Road
Watford WD17 1JJ
01923 237484
www.precisioncoping.co.uk
info@precisioncoping.co.uk



Laboratory:

Bin No:

Job No:

This custom-made appliance is for the use of:

Patient:

Shade:

Date Required:

Express Service:
Yes/No

Your Ref:

Lava

Procera

Other

18	17	16	15	14	13	12	11	21	22	23	24	25	26	27	28
48	47	46	45	44	43	42	41	31	32	33	34	35	36	37	38

Instructions:

N.B. For Lava the minimum coping thickness is 0.3mm, please consider the stump shades when requesting this option due to the translucency of Zirconia.

Enclosures:

Approved for
manufacture:

Date approved:

Date confirmed:

Technician:

Approved By:

This is a custom-made device intended for exclusive use by this patient. This device conforms to the relevant essential requirements set out in Annex I of the Medical Devices Directive. This statement does not apply to medical devices that have been repaired and/or refurbished for an individual patient's use.